

Phyco-Bubble

Airlift Photobioreactor





Varicon Aqua is the manufacturer of the **Phyco™** range of photobioreactors with more than 30 years' experience designing, constructing and deploying **algal photobioreactors**, high rate algal ponds and aquaculture algae production systems worldwide. As a global leader in this sector our close client engagement and supply chain relationships allow us to provide bespoke photobioreactor and cultivation systems for a broad range of applications. To date we have deployed over 270 photobioreactor systems across the globe, ranging in scale from 5 L to 400,000 L. Prominent products in the Phyco™ line include **Phyco-Flow™**, **Phyco-Pond** systems, the **Phyco-Lift** and **Phyco-Bubble**.

The **Phyco-Bubble** photobioreactor is an airlift system that provides a fully controllable laboratory scale column reactor system to the microalgae industry. The **Bubble** configuration is available in a range of incremental sizes from 5 - 15 L, and can be illuminated with LED or fluorescent Gro-Lux back lighting that can be tailored to specific photosynthetic requirements. The system is mixed by aeration which significantly reduces mechanical shear, and display excellent levels of turbulence and mass transfer. Standard fabrication is with fittings and connectors made from machined plastic, aluminium or stainless steel. The photostage is made from **Schott Duran® grade borosilicate**, producing a highly productive system with a superior lifespan.

Phyco-Bubble outline and options

All Phyco-Bubble™s can be fitted with a variety of ports for inlet/outlet and sensor functions dependent on requirements upon request.



Specifications and performance

Physical Considerations	
Volume/Height	5-15 L (12 L standard)
Tube OD	110 mm
Tube length	0.5-2 m (1.5 m standard)
Shape/mixing	Single tubular, mounted/aeration
Electricals	Lights, aerator, VAS control unit
Example Performance (details of other species available)	
Algae	<i>Chlorella sorokiniana</i> (species)
Density	1-2 g L ⁻¹
Productivity	0.1-0.3 g L ⁻¹ d ⁻¹
Max yield	1-2.5 g L ⁻¹
dO ₂	≤ 150 %
Engineering Performance (in house data)	
Liquid	0.1-0.2 m s ⁻¹ (max - velocity)
Reynolds	≥ 10 k (number)
Maximum	2,000 μ E m ⁻² s ⁻¹ (light)
Temperature	10-40 °C
Mass	~0.007 s ⁻¹ (transfer coefficient)

Phyco-Bubble configurations



Algae Production

In addition to our photobioreactor systems we also specialise in providing the auxiliary equipment for process monitoring and control, LED solutions and downstream processing. As a one-stop shop for algae production we have a full range of nutrient mixes from our **Cell-Hi All-in-One** line.

Our guarantee.

As with all Varicon Aqua systems we pride ourselves on the high level of after sales support we provide to our clients whether it's a systems issue or a technical production question. We are highly approachable. Come talk to us about your needs!

Contact: variconaqu@variconaqu.com - www.variconaqu.com

MAKING ALGAE EASY